

TaqMan Array Plates

Ready-to-use plates with preloaded TaqMan Assays



Introduction

Manual transfer of reagents into plates during routine experimental setups can be time-consuming and tedious. Daily preparation of such plates can negatively impact overall productivity. To ease the stress and minimize the error rate, Applied Biosystems™ TaqMan™ Array Plates are designed to reduce redundant pipetting steps and minimize loading errors, making them an ideal solution for fast and easy experimental setup.

TaqMan Array Plates contain high-quality Applied Biosystems™ TaqMan™ Gene Expression Assays (TaqMan™ probe and PCR primer sets) in a convenient, dried-down, 96- or 384-well plate format. Extensive tests have found no difference in performance between TaqMan Assays supplied in solution in single tubes and assays supplied dried in a TaqMan Array Plate (Figure 1).

- **Convenient setup**—dried TaqMan Assays in a 96- or 384-well plate; just add master mix and your cDNA sample
- **Flexible format**—choose from preconfigured to customizable TaqMan Array Plates in fast or standard formats
- **Easy data analysis**—Applied Biosystems™ analysis modules help provide rapid cloud-based data interpretation




Figure 1. Linear dynamic range of TaqMan Array Plates (standard) vs. TaqMan Assays.

TaqMan Array Plates have the same performance as standard TaqMan Assays (wet) on a 96-well plate. The linear dynamic range is maintained, even down to low cDNA inputs, showing the assays are unaffected by the process of drying down the TaqMan Array Plate. Each plate contained one assay (18S Hs99999901_s1) and eight replicates for each cDNA input.

TaqMan Array Plates include enough TaqMan Assay in each well for a 20 µL reaction (96-well standard plate), 10 µL reaction (96-well fast plate), or 5–10 µL reaction (384-well plate). These preloaded plates can help accelerate research in a variety of application areas, including microarray validation and pathway studies (Table 1).

Table 1. Examples of published research using TaqMan Array Plates for a variety of applications.

Research focus	Title	Publication
Cancer	Benzene-poly-carboxylic acid complex, a novel anti-cancer agent induces apoptosis in human breast cancer cells	<i>PLoS One</i> 9:e85156 (2014)
Autophagy	L450W and Q455K Col8a2 knock-in mouse models of Fuchs endothelial corneal dystrophy show distinct phenotypes and evidence for altered autophagy	<i>Invest Ophthalmol Vis Sci</i> 54:1887 (2013)
Infectious diseases	Distinct and overlapping genomic profiles and antiviral effects of Interferon-λ and -α on HCV-infected and non-infected hepatoma cells	<i>J Viral Hepat</i> 19:843 (2012)
Apoptosis	Ribonuclease binase apoptotic signature in leukemic Kasumi-1 cells	<i>Biochimie</i> 95:1344 (2013)
Angiogenesis	Modulation of angiogenesis by genetic manipulation of ATF4 in mouse model of oxygen-induced retinopathy	<i>Invest Ophthalmol Vis Sci</i> 54:5995 (2013)

Versatile formats

TaqMan Array Plates are available in various formats to meet your laboratory’s needs (Table 2).

- **Preconfigured (fixed-content) TaqMan Array Plates**—fixed content of the most sought-after predefined gene panels, categorized by specific disease, pathway, or biological process (Table 3)
- **Flexible-content TaqMan Array Plates**—select preconfigured pathway panels and modify with predesigned TaqMan Assays to suit your needs
- **Custom TaqMan Array Plates**—can be readily designed from more than 1.3 million predesigned TaqMan Gene Expression Assays available for 24 species, including some pathogens. Ordering is easy with the online custom TaqMan Array Plate configuration tool, which helps you find and select genes and assays. Custom TaqMan Array Plates are available in 9 different formats (Table 4).

- **Custom formatting service**—design your own TaqMan Array Plate configured to your specifications (e.g., Applied Biosystems™ TaqMan™ SNP Genotyping Assays, Applied Biosystems™ TaqMan™ Copy Number Assays, or combined mRNA/miRNA assays). Custom plating services are available for 96-well standard, 96-well fast, and 384-well plates.

Easy workflow, superior sensitivity

Save a significant amount of time setting up your experiments—just add master mix and cDNA sample, and begin cycling. The standard plate format provides real-time PCR results in about 1.5 hours, whereas the fast plate format typically delivers results in only 30 to 45 minutes.

Furthermore, you can use the optional Applied Biosystems™ TaqMan™ PreAmp Master Mix and Applied Biosystems™

Table 2. TaqMan Array Plate formats and application compatibility.

	Preconfigured TaqMan Array Plates	Flexible-content TaqMan Array Plates	Custom TaqMan Array Plates	Custom formatting service
Definition	Fixed, predefined gene panels	Modifiable predefined gene panels	Configurable predesigned assays	Client-requested designs
Web page	thermofisher.com/taqmanarrays	thermofisher.com/flexiblepanels	thermofisher.com/arrayplates	thermofisher.com/customformattingsservice
Application supported				
Gene expression	√ (96-well)	√ (96-well)	√ (96-well)	√ (96-, 384-well)
MicroRNA				√ (96-, 384-well)
Genotyping				√ (96-, 384-well)

Applied Biosystems™ qPCR analysis modules are free, easy-to-use cloud-based data analysis tools with analysis module options for comparative C_t analysis, also known as relative quantification (RQ), and standard curve analysis. It provides integrated analysis of multiple data sets, while offering new functionalities such as an online file storage system, flexible plate setup, analysis groups, and robust visualization to place your data fully in your control (Figure 2).



Figure 2. Example of data visualization with Applied Biosystems qPCR analysis modules. (A) Relative quantification (RQ) plot view. **(B)** Heat map view for an RQ study.

Table 3. Preconfigured TaqMan Array Plates. This is a partial list; catalog numbers in this table are for human unless otherwise noted. For the complete list, visit thermofisher.com/tagmanarrays

Preconfigured (fixed-content) TaqMan Array Plates	Cat. No. (Standard 96-well)	Cat. No. (Fast 96-well)
Antioxidant mechanisms	4414119	4418764
Apoptosis	4414072	4418717
Drug transporters	4414118	4418763
Endogenous controls	4396840	4426700
Extracellular matrix and adhesion molecules	4414133	4418778
Inflammation	4414074	4418719
Lipid-regulated genes	4415461*	4418846*
Molecular mechanisms of cancer	4414161	4418806
Phagocytosis of microbes	4414178	4418823
Rheumatoid arthritis pathogenesis	4414170	4418815
Ubiquitin-proteasome proteolysis	4414198	4418843
WNT pathway	4414100	4418745

*Cat. No. is for mouse.

Table 4. Custom TaqMan Array Plate format options for 96-well plates.

Format	No. of assays + controls	No. of samples	Minimum order	Cat. No. (standard)	Cat. No. (fast)
TaqMan Array Plate 8	7 + 1	12	6	4413266	4413263
TaqMan Array Plate 16	15 + 1	6	6	4413264	4413261
TaqMan Array Plate 16 Plus	12 + 4	6	6	4413265	4413262
TaqMan Array Plate 32	31 + 1	3	6	4391528	4413259
TaqMan Array Plate 32 Plus	28 + 4	3	6	4391529	4413260
TaqMan Array Plate 48	47 + 1	2	6	4391526	4413257
TaqMan Array Plate 48 Plus	44 + 4	2	6	4391527	4413258
TaqMan Array Plate 96	95 + 1	1	6	4391524	4413255
TaqMan Array Plate 96 Plus	92 + 4	1	6	4391525	4413256

Ordering information—related products

Product	Quantity	Cat. No.
Preamplification reagents		
TaqMan PreAmp Master Mix (2X)	40 reactions	4391128
TaqMan PreAmp Master Mix (2X)	200 reactions	4488593
Custom TaqMan PreAmp Pools	250 reactions	4441856

For compatible TaqMan Master Mixes, visit thermofisher.com/taqmanmastermixes

For 384-well plate options, visit thermofisher.com/customformatting-service

Find out more at thermofisher.com/taqmanarrays

ThermoFisher
SCIENTIFIC

TaqMan Array Cards

Load 384 reactions accurately in 10 minutes without robotics


- **Easy and rapid setup**—load 384 wells in approximately 10 minutes, without expensive liquid-handling robotics
- **Highly reproducible results**—ideal for low-expressing genes or precious samples
- **Flexible format designs**—choose from preconfigured panels or a full custom design
- **Easy data analysis**—enables rapid and accurate analysis across a large number of genes and samples

Introduction

As experimental throughput in research laboratories (both clinical and nonclinical) increases, so does the need for simple and efficient ways to run tests. Researchers are often left to choose between cumbersome, time-consuming manual approaches and costly automation. Applied Biosystems™ TaqMan™ Array Cards provide an alternative that enables you to achieve highly reproducible and sensitive results with higher throughput, but without the expense of liquid-handling robotics.

Customized for your application

Widely cited in publications (Table 1), TaqMan Array Cards are ideal for medium-throughput validation and screening studies, whether your research involves cancer, stem cells, inflammation, or infectious diseases. Use the cards for validating tens or hundreds of initial hits generated from microarrays or next-generation sequencing, or use TaqMan Array Cards for screening biomarkers and toxicology panels, or for analyzing pathways, target classes, and complete disease sets.



A TaqMan Array Card is a 384-well microfluidic card designed for performing 384 simultaneous real-time PCR reactions without the need for expensive liquid-handling automation. TaqMan Array Cards are preloaded with dried-down high-quality TaqMan Assays (TaqMan probes and PCR primer sets), ready for 1 to 8 samples to be run in parallel against 12 to 384 Applied Biosystems™ TaqMan™ Gene Expression Assay targets (including a manufacturing control) (Figure 1). Because of their design, TaqMan Array Cards make it easy to produce consistent results with low variability across multiple users and laboratories. You can rely on TaqMan Array Cards to help you quickly achieve highly reproducible and sensitive results.

Table 1. Examples of published research using TaqMan Array Cards for a variety of applications.

Research focus	Title	Publication
Stem cells	Different isolation methods alter the gene expression profiling of adipose-derived stem cells	<i>Int J Med Sci</i> 11:391 (2014)
	Stem cells expanded from the human embryonic hindbrain stably retain regional specification and high neurogenic potency	<i>J Neurosci</i> 33:12407 (2013)
	Isolation, characterization, and gene expression analysis of Wharton's jelly-derived mesenchymal stem cells under xeno-free culture conditions	<i>Stem Cells Cloning</i> 4:39 (2011)
Cancer	Differentially expressed miRNAs in Ewing sarcoma compared to mesenchymal stem cells: low miR-31 expression with effects on proliferation and invasion	<i>PLoS One</i> 9:e93067 (2014)
	Oncogenic micro-RNAs and renal cell carcinoma	<i>Front Oncol</i> 4:49 (2014)
	Gene expression analysis of so called Asian dust extracts in human acute myeloid leukemia cells	<i>Toxicol Res</i> 26:21 (2010)
Inflammation	Expression of genes related to anti-inflammatory pathways are modified among farmers' children	<i>PLoS One</i> 9:e91097 (2014)
	Resolution of central nervous system astrocytic and endothelial sources of <i>CCL2</i> gene expression during evolving neuroinflammation	<i>Fluids Barriers CNS</i> 11:6 (2014)
	Inflammatory activation is associated with a reduced glucocorticoid receptor alpha/beta expression ratio in monocytes of inpatients with melancholic major depressive disorder	<i>Transl Psychiatry</i> 4:e344 (2014)
Toxicity	Selenoprotein P regulates 1-(4-Chlorophenyl)-benzo-2,5-quinone-induced oxidative stress and toxicity in human keratinocytes	<i>Free Radic Biol Med</i> 65:70 (2013)
	Expression profiling of selected genes of toxication and detoxication pathways in peripheral blood lymphocytes as a biomarker for predicting toxicity of environmental chemicals	<i>Int J Hyg Environ Health</i> 216:645 (2013)
Infectious diseases	Detection and characterization of mycoplasma pneumoniae during an outbreak of respiratory illness at a university	<i>J Clin Microbiol</i> 52:849 (2014)
	Optimization of multiple pathogen detection using the TaqMan Array Card: application for a population-based study of neonatal infection	<i>PLoS One</i> 8:e66183 (2013)
	Field evaluation of TaqMan Array Card (TAC) for the simultaneous detection of multiple respiratory viruses in children with acute respiratory infection	<i>J Clin Virol</i> 57:254 (2013)




Figure 1. TaqMan Array Card and custom formats. (A) The TaqMan Array Card, designed to be used with the Applied Biosystems™ QuantStudio™ 7, QuantStudio™ 12 Flex, ViiA™ 7, or 7900HT Real-Time PCR System, is a real-time PCR “lab on a chip” containing 384 wells connected by a series of microfluidic channels. Each loading port is connected to 48 microwells of ~1.5 µL containing selected dried TaqMan Assays. **(B)** Each port can be loaded with the same or different samples, allowing analysis of 1 to 8 samples per card. Custom TaqMan Array Cards are available in 10 different formats with 12, 16, 24, 32, 48, 64, 96 (2 choices), 192, and 384 assays.

Simple and straightforward loading

Using the TaqMan Array Card is fast and simple (Figure 2). The card has 8 sample-loading ports, each connected to a set of 48 reaction wells. Simply pipette your cDNA sample premixed with Applied Biosystems™ TaqMan™ Universal Master Mix II into each port, briefly centrifuge to fill the wells, and seal the card to close the wells.* Within 10 minutes, your card is ready to run on the Applied Biosystems QuantStudio 7, QuantStudio 12K Flex, ViiA 7, or 7900HT Real-Time PCR System equipped with an Applied Biosystems™ TaqMan™ Array block. This streamlines reaction setup, saves time, and reduces labor-intensive pipetting steps. No need for hours of training to master the preparation of the card—a simple 5-minute video provides step-by-step instructions to make learning easy.




Figure 2. TaqMan Array Card workflow.*

Versatile formats

TaqMan Array Cards are available in various formats to meet your laboratory's needs (Table 2).

- **Preconfigured (fixed-content) TaqMan Array Cards**—choose premade cards that contain the most widely used predefined gene expression panels (categorized by specific disease, pathway, or biological process) or miRNA panels (Table 3)
- **Flexible-content TaqMan Array Cards**—select pre-configured pathway panels and modify assay content to suit your needs
- **Custom TaqMan Array Cards**—design them yourself, choosing from any of over 1.3 million predesigned Applied Biosystems™ TaqMan™ Gene Expression Assays available for 24 species, or over 3,000 Applied Biosystems™ TaqMan™ miRNA assays. Ordering is easy with the online Applied Biosystems™ Custom TaqMan Array Card configuration tool, which helps you find and select genes and assays. Custom TaqMan Array Cards are available in 10 different formats (Figure 1).
- **Custom formatting service**—design your own TaqMan Array Card configured to your specifications (e.g., Applied Biosystems™ TaqMan™ SNP Genotyping Assays, Applied Biosystems™ TaqMan™ Copy Number Assays, combined

Table 2. TaqMan Array Card formats and application compatibility.

	Preconfigured TaqMan Array Cards	Flexible-content TaqMan Array Cards	Custom TaqMan Array Cards	Custom formatting service
Definition	Fixed, predefined gene expression or miRNA panels (Table 3)	Modifiable, predefined panels	Configurable predesigned assays	Customized designs
Web page	thermofisher.com/taqmanarrays	thermofisher.com/flexiblepanels	thermofisher.com/arraycards	thermofisher.com/customformatting-service
Application supported				
Gene expression	✓	✓	✓	✓
MicroRNA	✓		✓	✓
Genotyping				✓

*TaqMan Array Card centrifuge compatibility chart available at thermofisher.com/centrifugecompatibility

Table 3. Preconfigured (fixed-content) TaqMan Array Cards.

Gene expression panels			
Panel	Cat. No. (Human)	Cat. No. (Mouse)	Cat. No. (Rat)
ABC transporter	4378700		
Alzheimer's	4378713	4378714	
Angiogenesis	4378710		
Apoptosis	4378701		
Endogenous control	4367563	4378702	4378704
GPCR	4367785	4378703	4378709
Hepatocarcinogenicity			4465484
Immune	4370573	4367786	
Inflammation	4378707		4378708
Nuclear receptor	4379961		
Phosphodiesterase	4378705		4378706
Protein kinase	4367784		
Stem cell pluripotency	4385344	4385363	
Viral pathogen screening	4465485		

miRNA panels		
Panel	Cat. No. (Human)	Cat. No. (Rodent)
TaqMan Array MicroRNA A+B Cards Set v3.0	4444913	4444909
TaqMan Array MicroRNA A Cards v2.0	4398965	4398967
TaqMan Array MicroRNA B Cards v3.0	4444910	4444899

assays for mRNA and miRNA or other species)

Powerful data analysis

Applied Biosystems™ qPCR Analysis Modules are free, easy-to-use data analysis tools for comparative C_t analysis, also known as relative quantification (RQ), and standard curve analysis. It provides integrated analysis of multiple data sets, while offering new functionalities such as an online file storage system, flexible plate setup, analysis groups, and robust visualization to place your data fully in your control (Figure 3). Analyze up to 500 TaqMan Array Cards in one study with Applied Biosystems analysis modules.




Figure 3. Example of data visualization with Applied Biosystems analysis modules. (A) Relative quantification (RQ) plot view. (B) Heat map view for an RQ study.

High-quality results

Results from TaqMan Array Cards are highly reproducible, both within and across individual array cards. Uniform distribution of samples to the well chambers with minimal handling helps to significantly reduce the risk of introducing variability or cross-contamination. Rely on the high precision and sensitivity of the results that TaqMan Array Cards can help you achieve, especially for precious samples with minimally expressed targets (Figure 4). You may also use the optional Applied Biosystems™ TaqMan™ PreAmp Master Mix and Custom PreAmp Pools to generate a comprehensive expression profile with a small sample input—as little as 1 ng of total RNA. Preamplification can enhance the ability to detect low-abundance RNA targets and help your precious sample last for many more real-time PCR runs.

TaqMan Array Card specifications

Loading time	<10 minutes
Volume per well	~1.5 µL
Nucleic acid template	30–1,000 ng
Loading volume	800 µL/card (100 µL per port, cDNA and master mix combined)
Run mode support	Standard and fast modes
Assay throughput	11–380 assays/card
Sample throughput	1–8 samples/card

A

Dilution	Concentration pg/well	Card 1		Card 2		Card 3		Cards 1-3	
		C _t mean	C _t SD	C _t mean	C _t SD	C _t mean	C _t SD	C _t mean	C _t SD
1	10,000	6.32	0.13	6.57	0.26	6.39	0.11	6.43	0.13
2	1,000	9.86	0.24	9.81	0.07	9.71	0.15	9.79	0.08
3	100	13.23	0.08	13.17	0.16	13.15	0.04	13.18	0.04
4	10	16.58	0.17	16.58	0.09	16.45	0.20	16.54	0.08
5	1	20.00	0.13	20.02	0.28	19.92	0.24	19.98	0.05
6	0.1	23.37	0.13	23.34	0.14	23.29	0.17	23.33	0.04
7	0.01	26.80	0.14	26.99	0.15	26.80	0.17	26.86	0.11
8	0.001	30.51	0.20	30.44	0.33	30.36	0.57	30.44	0.07




Figure 4. Achieve high reproducibility with broad dynamic range with TaqMan Array Cards. (A) Amplification of the 18S gene using 1 µg–0.1 pg of cDNA per well was performed on 3 different TaqMan Array Cards to evaluate card-to-card reproducibility. The table shows average C_t values and standard deviations (SD) for each dilution for each card. Card-to-card standard deviation of the C_t mean is 0.13 or less for all dilutions, showing good reproducibility at both low and high target concentrations. **(B)** The amplification plot for card 1 (n = 3).

Ordering information

Product	Cat. No.	
Custom TaqMan Array Cards		
Custom Gene Expression TaqMan Array Card - Format 12	4342247	
Custom Gene Expression TaqMan Array Card - Format 16	4346798	
Custom Gene Expression TaqMan Array Card - Format 24	4342249	
Custom Gene Expression TaqMan Array Card - Format 32	4346799	
Custom Gene Expression TaqMan Array Card - Format 48	4342253	
Custom Gene Expression TaqMan Array Card - Format 64	4346800	
Custom Gene Expression TaqMan Array Card - Format 96a	4342259	
Custom Gene Expression TaqMan Array Card - Format 96b	4342261	
Custom Gene Expression TaqMan Array Card - Format 192	4346802	
Custom Gene Expression TaqMan Array Card - Format 384	4342265	
Custom TaqMan Array MicroRNA Cards - Format 12	4449135	
Custom TaqMan Array MicroRNA Cards - Format 16	4449136	
Custom TaqMan Array MicroRNA Cards - Format 24	4449137	
Custom TaqMan Array MicroRNA Cards - Format 32	4449138	
Custom TaqMan Array MicroRNA Cards - Format 48	4449139	
Custom TaqMan Array MicroRNA Cards - Format 64	4449140	
Custom TaqMan Array MicroRNA Cards - Format 96a	4449141	
Custom TaqMan Array MicroRNA Cards - Format 96b	4449142	
Custom TaqMan Array MicroRNA Cards - Format 192	4449143	
Custom TaqMan Array MicroRNA Cards - Format 384	4449144	
PreAmp Pool Reagents	Quantity	
TaqMan PreAmp Master Mix (2X)	40 reactions	4391128
TaqMan PreAmp Master Mix (2X)	200 reactions	4488593
Custom TaqMan PreAmp Pools	250 reactions	4441856
TaqMan Array Card Blocks†		
QuantStudio 12K Flex TaqMan Array Card Block		4453546
QuantStudio 7 TaqMan Array Card Block		4453546
ViiA 7 TaqMan Array Card Block		4453546
7900HT TaqMan Array Card Block		4329012

†Applied Biosystems™ TaqMan™ Array Card Instrument Blocks include a sample block, TaqMan Array Card sealer, custom centrifuge buckets with adaptors, getting started guide, and a chemistry installation kit.

To find compatible TaqMan master mixes, go to thermofisher.com/taqmanmm

To find out more, go to thermofisher.com/taqmanarrays

ThermoFisher
SCIENTIFIC